

LIN28 antibody



Catalog Number: 107377

Product name

LIN28 antibody

Immunogen

Specificity

Human; other species not tested.

Antibody description

LIN28 Mouse Monoclonal antibody. Positive WB detected in NCCIT cell. Observed molecular weight by Western-blot: 28 kDa

Preparation

This antibody was obtained by immunization of LIN28 recombinant protein (Accession Number: NM_024674). Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



NCCIT cell were subjected to SDS PAGE followed by western blot with (LIN28 Antibody) at dilution of 1:1000