

IL37 antibody

Catalog Number: 107400

General Information

Product name:

IL37 antibody

Immunogen:

Applications:

ELISA, WB

Specificity:

Human; other species not tested.

Antibody description:

IL37 Mouse Monoclonal antibody. Positive WB detected in Transfected HEK-293 cells.

Properties

Preparation:

This antibody was obtained by immunization of IL37 recombinant protein (Accession Number: BC020637). Purification method: Protein G purified.

Formulation:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage:

Store at -20°C. DO NOT ALIQUOT

Clonality:

Monoclonal


Ig Type: Mouse IgG2a

Applications

Recommended Dilution:

WB: 1:1000-1:10000

Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:107400(IL1F7 Antibody) at dilution of 1:2000