

# MAPK13 antibody



Catalog Number: 107415

## Product name

Monoclonal

MAPK13 antibody

## Ig Type

## Specificity

Mouse IgG1

Human, Mouse; other species not tested.

## Applications

## Antibody description

ELISA, WB

MAPK13 Mouse Monoclonal antibody. Positive WB detected in HepG2 cells, HeLa cells, MCF-7 cells. Observed molecular weight by Western-blot: 42kd

## Dilutions

Recommended Dilution:

## Preparation


WB: 1:500-1:5000

This antibody was obtained by immunization of MAPK13 recombinant protein (Accession Number: NM\_002754). Purification method: Protein G purified.

## Validations

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



## Storage

Store at -20°C. DO NOT ALIQUOT

HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:107415(MAPK13 antibody) at dilution of 1:1000

## Clonality