

MCRS1 antibody



Catalog Number: 107421

Product name

MCRS1 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

MCRS1 Mouse Monoclonal antibody. Positive WB detected in Skov-3 cells. Observed molecular weight by Western-blot: 55+110KD

Preparation

This antibody was obtained by immunization of MCRS1 recombinant protein (Accession Number: NM_006337). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2a

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



SKOV-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:107421(MCRS1 antibody) at dilution of 1:1000