

p38 MAPK antibody



Catalog Number: 107443

Product name

p38 MAPK antibody

Specificity

Human, Pig; other species not tested.

Antibody description

p38 MAPK Mouse Monoclonal antibody. Positive WB detected in human heart tissue, pig heart tissue. Observed molecular weight by Western-blot: 38-42 kDa

Preparation

This antibody was obtained by immunization of p38 MAPK recombinant protein (Accession Number: NM_139012). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2b

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:107443(p38 MAPK Antibody) at dilution of 1:2000