

# p38 MAPK antibody



Catalog Number: 107443

## Product name

p38 MAPK antibody

## Immunogen

[Human p38 MAPK Recombinant protein \(His tag\)](#)

## Specificity

Human, Pig; other species not tested.

## Antibody description

p38 MAPK Mouse Monoclonal antibody. Positive WB detected in human heart tissue, pig heart tissue. Observed molecular weight by Western-blot: 38-42 kDa

## Preparation

This antibody was obtained by immunization of p38 MAPK recombinant protein (Accession Number: NM\_139012). Purification method: Protein A purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2b


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:1000-1:10000

## Validations



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:107443(p38 MAPK Antibody) at dilution of 1:2000