

RXRG antibody



Catalog Number: 107515

Product name

RXRG antibody

Immunogen

[Human RXRG Recombinant protein \(GST tag\)](#)

Specificity

Human; other species not tested.

Antibody description

RXRG Mouse Monoclonal antibody. Positive WB detected in A549 cells. Observed molecular weight by Western-blot: 52 kDa

Preparation

This antibody was obtained by immunization of RXRG recombinant protein (Accession Number: NM_006917). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2b


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:1000-1:10000

Validations



A549 cells were subjected to SDS PAGE followed by western blot with Catalog No:107515(RXRG Antibody) at dilution of 1:2000