

# RXRG antibody



Catalog Number: 107515

## Product name

RXRG antibody

## Immunogen

[Human RXRG Recombinant protein \(GST tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

RXRG Mouse Monoclonal antibody. Positive WB detected in A549 cells. Observed molecular weight by Western-blot: 52 kDa

## Preparation

This antibody was obtained by immunization of RXRG recombinant protein (Accession Number: NM\_006917). Purification method: Protein A purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2b


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:1000-1:10000

## Validations



A549 cells were subjected to SDS PAGE followed by western blot with Catalog No:107515(RXRG Antibody) at dilution of 1:2000