

WTAP antibody



Catalog Number: 107544

Product name

WTAP antibody

WB: 1:1000-1:10000

IHC: 1:20-1:200

Specificity

Human, Mouse; other species not tested.

IF: 1:20-1:200

Antibody description

WTAP Mouse Monoclonal antibody. Positive WB detected in HEK-293 cells, HeLa cells, HepG2 cells, human testis tissue, Jurkat cells, mouse testis tissue. Positive IF detected in A431 cells. Positive IHC detected in human cervical cancer tissue. Observed molecular weight by Western-blot: 50 kDa

Preparation

This antibody was obtained by immunization of WTAP recombinant protein (Accession Number: NM_152858). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1


Applications

ELISA, IHC, IF, WB


Dilutions

Recommended Dilution:


Validations



Immunohistochemical of paraffin-embedded human cervical cancer using Catalog No:107544(WTAP antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human cervical cancer using Catalog No:107544(WTAP antibody) at dilution of 1:50




(under 40x lens)

WTAP antibody



Catalog Number: 107544

HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:107544(WTAP Antibody) at dilution of 1:2000



Immunofluorescent analysis of A431 cells using Catalog No:107544(WTAP Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Mouse IgG(H+L)