

SLAM/CD150 antibody



Catalog Number: 107588

Product name

Monoclonal

SLAM/CD150 antibody

Ig Type

Specificity

IgM

Human; other species not tested.

Applications

Antibody description

ELISA, WB

SLAM/CD150 mouse Monoclonal antibody.
Positive WB detected in Raji cells, HUVEC cells,
Jurkat cells. Observed molecular weight by
Western-blot: 80 kDa

Dilutions


Recommended Dilution:

WB: 1:1000-1:10000

Preparation

This antibody was obtained by immunization of
SLAM/CD150 recombinant protein (Accession
Number: BC012602). Purification method:
Caprylic acid/ammonium sulfate precipitation.

Validations



Formulation

PBS with 0.1% sodium azide and 50% glycerol pH
7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Raji cells were subjected to SDS PAGE followed by
western blot with Catalog No:107588(SLAMF1
Antibody) at dilution of 1:2000