

SULT4A1 antibody



Catalog Number: 107613

Product name

Monoclonal

SULT4A1 antibody

Ig Type

Specificity

Mouse IgG1

Human; other species not tested.

Applications

Antibody description

ELISA, WB

SULT4A1 Mouse Monoclonal antibody. Positive WB detected in Human brain. Observed molecular weight by Western-blot: 33kDa

Dilutions


Recommended Dilution:

WB: 1:500-1:5000

Preparation

This antibody was obtained by immunization of SULT4A1 recombinant protein (Accession Number: NM_014351). Purification method: Protein G purified.

Validations



Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:107613(SULT4A1 antibody) at dilution of 1:1000