

TH antibody



Catalog Number: 107622

Product name

TH antibody

Specificity

Human; other species not tested.

Antibody description

TH Mouse Monoclonal antibody. Positive WB detected in SH-SY5Y cells, PC-12 cells. Observed molecular weight by Western-blot: 55 kDa

Preparation

This antibody was obtained by immunization of TH recombinant protein (Accession Number: NM_000360). Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:107622(TH Antibody) at dilution of 1:2000