

ABI2 antibody



Catalog Number: 107817

Product name

ABI2 antibody

WB: 1:500-1:5000

IHC: 1:20-1:200

Specificity


Human, Mouse, Rat; other species not tested.

IF: 1:10-1:100

Validations

Antibody description

ABI2 Rabbit Polyclonal antibody. Positive IF detected in Hela cells. Positive IHC detected in human brain tissue, human testis tissue. Positive WB detected in mouse heart tissue, HEK-293 cells, mouse brain tissue, mouse testis tissue, mouse uterus tissue. Observed molecular weight by Western-blot: 60-70 kDa



Preparation

This antibody was obtained by immunization of ABI2 recombinant protein (Accession Number: NM_005759). Purification method: Antigen affinity purified.

mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:107817(ABI2 antibody) at dilution of 1:400

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications


ELISA, WB, IF, IHC

Dilutions

Recommended Dilution:




Immunofluorescent analysis of Hela cells, using ABI2 antibody Catalog No:107817 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).



ABI2 antibody

Catalog Number: 107817

Immunohistochemistry of paraffin-embedded human brain tissue slide using Catalog No:107817(ABI2 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human brain tissue slide using Catalog No:107817(ABI2 Antibody) at dilution of 1:50 (under 40x lens)