

ACAA1 antibody

Catalog Number: 107823

Product name

ACAA1 antibody

WB: 1:500-1:5000

IHC: 1:20-1:200

Specificity


Human, Mouse, Rat; other species not tested.

IF: 1:10-1:100

Antibody description

ACAA1 Rabbit Polyclonal antibody. Positive WB detected in L02 cells, A549 cells, human lung tissue. Positive IF detected in HeLa cells. Positive IHC detected in human thyroid cancer tissue, human liver cancer tissue, human prostate cancer tissue. Observed molecular weight by Western-blot: 41 kDa

Validations




Preparation

This antibody was obtained by immunization of ACAA1 recombinant protein (Accession Number: NM_001607). Purification method: Antigen affinity purified.

L02 cells were subjected to SDS PAGE followed by western blot with Catalog No:107823(ACAA1 antibody) at dilution of 1:1000

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Immunohistochemistry of paraffin-embedded human thyroid cancer tissue slide using Catalog No:107823(ACAA1 Antibody) at dilution of 1:50 (under 10x lens)

Ig Type


Rabbit IgG

Applications

ELISA, WB, IF, IHC

Dilutions


Recommended Dilution:



ACAA1 antibody

Catalog Number: 107823

Immunohistochemistry of paraffin-embedded human thyroid cancer tissue slide using Catalog No:107823(ACAA1 Antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of HeLa cells using Catalog No:107823(ACAA1 Antibody) at dilution of 1:25 and Rhodamine-Goat anti-Rabbit IgG