

# AGPAT3 antibody



Catalog Number: 107832

## Product name

AGPAT3 antibody

## Immunogen

[Human AGPAT3 Recombinant protein \(His tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

AGPAT3 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, rat testis tissue. Observed molecular weight by Western-blot: 36-40 kDa

## Preparation

This antibody was obtained by immunization of AGPAT3 recombinant protein (Accession Number: XM\_017028410). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:500-1:5000

## Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:107832(AGPAT3 Antibody) at dilution of 1:1000