

AMBRA1 antibody



Catalog Number: 107988

Product name

AMBRA1 antibody

WB: 1:200-1:1000

IHC: 1:20-1:200

Specificity


Human, Mouse, Rat; other species not tested.

IF: 1:10-1:100

Antibody description

AMBRA1 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, HepG2 cells, HT-1080 cells, human brain tissue, MCF7 cells. Positive IHC detected in human liver cancer tissue. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 130-135kd

Validations




Preparation

This antibody was obtained by immunization of AMBRA1 recombinant protein (Accession Number: XM_005253014). Purification method: Antigen affinity purified.

mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:107988(AMBRA1 antibody) at dilution of 1:300

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Immunohistochemical of paraffin-embedded human liver cancer using Catalog No:107988(AMBRA1 antibody) at dilution of 1:50 (under 10x lens)

Ig Type


Rabbit IgG

Applications

ELISA, IHC, WB, IF

Dilutions

Recommended Dilution:



AMBRA1 antibody



Catalog Number: 107988

Immunofluorescent analysis of HepG2 cells using Catalog No:107988(AMBRA1 Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)