

# Angiopoietin 1 antibody



Catalog Number: 108051

## Product name

Angiopoietin 1 antibody

## Immunogen

[Human ANGPT1 Recombinant protein \(His tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

Angiopoietin 1 Rabbit Polyclonal antibody. Positive IF detected in HeLa cells. Positive WB detected in mouse heart tissue, HepG2 cells, mouse liver tissue, rat heart tissue. Observed molecular weight by Western-blot: 70 kDa

## Preparation

This antibody was obtained by immunization of Angiopoietin 1 recombinant protein (Accession Number: NM\_001146). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

IF: 1:10-1:100

## Validations



mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:108051(ANGPT1 antibody) at dilution of 1:1000

Immunofluorescent analysis of HeLa cells using Catalog No:108051(ANGPT1 Antibody) at dilution of 1:25 and Rhodamine-Goat anti-Rabbit IgG