

Annexin A2 antibody



Catalog Number: 108093

Product name

Annexin A2 antibody

Immunogen

[Human Annexin A2 Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Annexin A2 Rabbit Polyclonal antibody. Positive IHC detected in human liver tissue, human colon cancer tissue, human oesophagus cancer tissue. Positive IF detected in HeLa cells, NIH/3T3 cells. Positive IP detected in HeLa cells. Positive WB detected in HeLa cells, K-562 cells, mouse placenta tissue, NIH/3T3 cells, rat uterus tissue. Observed molecular weight by Western-blot: 39 kDa

Preparation

This antibody was obtained by immunization of Annexin A2 recombinant protein (Accession Number: XM_017022090). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IP, WB, IF, IHC

Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IP: 1:200-1:2000


IHC: 1:50-1:500

IF: 1:50-1:500

Validations




HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:108093(Annexin A2 Antibody) at dilution of 1:2000




Immunohistochemistry of paraffin-embedded human liver tissue slide using Catalog No:108093(Annexin A2 Antibody) at dilution of 1:200 (under 10x lens).

Annexin A2 antibody


Catalog Number: 108093



Immunohistochemistry of paraffin-embedded human liver tissue slide using Catalog No:108093(Annexin A2 Antibody) at dilution of 1:200 (under 40x lens).



Immunofluorescent analysis of (-20oc Ethanol) fixed HeLa cells using Catalog No:108093(Annexin A2 Antibody) at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)



IP Result of anti-ANXA2 (IP:Catalog No:108093, 3ug; Detection:Catalog No:108093 1:500) with HeLa cells lysate 3000ug.