

# ATPIF1 antibody

Catalog Number: 108324

## Product name

ATPIF1 antibody

IP: 1:200-1:2000

IF: 1:50-1:500

## Specificity

Human, Mouse; other species not tested.

## Validations

## Antibody description

ATPIF1 Rabbit Polyclonal antibody. Positive IF detected in HeLa cells. Positive IP detected in HeLa cells. Positive WB detected in HeLa cells. Observed molecular weight by Western-blot: 12 kDa




## Preparation

This antibody was obtained by immunization of ATPIF1 recombinant protein (Accession Number: NM\_016311). Purification method: Antigen affinity purified.

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:108324(ATPIF1 antibody) at dilution of 1:500

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

IP Result of anti-ATPIF1 (IP:Catalog No:108324, 4ug; Detection:Catalog No:108324 1:500) with HeLa cells lysate 3440ug.

## Ig Type

Rabbit IgG


## Applications

ELISA, WB, IP, IF

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000



Immunofluorescent analysis of (10% Formaldehyde) fixed HeLa cells using Catalog

# ATPIF1 antibody



Catalog Number: 108324

---

No:108324(ATPIF1 Antibody) at dilution of 1:50

and Alexa Fluor 488-conjugated AffiniPure Goat  
Anti-Rabbit IgG(H+L)