

# BCL6 antibody

Catalog Number: 108429

## Product name

BCL6 antibody

IP: 1:500-1:5000

IHC: 1:50-1:500


## Specificity

Human; other species not tested.

## Validations

## Antibody description

BCL6 Rabbit Polyclonal antibody. Positive WB detected in Raji cells. Positive IP detected in Raji cells. Positive IHC detected in human lymphoma tissue, human tonsillitis tissue. Observed molecular weight by Western-blot: 80 kDa




## Preparation

This antibody was obtained by immunization of BCL6 recombinant protein (Accession Number: NM\_001130845). Purification method: Antigen affinity purified.

Immunohistochemistry of paraffin-embedded human lymphoma tissue slide using Catalog No:117102(BCL6 Antibody) at dilution of 1:400 (under 10x lens). heat mediated antigen retrieved with Sodium Citrate buffer (pH6).

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP, IHC


## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Immunohistochemistry of paraffin-embedded human lymphoma tissue slide using Catalog No:117102(BCL6 Antibody) at dilution of 1:400 (under 40x lens). heat mediated antigen retrieved

with Sodium Citrate buffer (pH6).



# BCL6 antibody

Catalog Number: 108429

Raji cells were subjected to SDS PAGE followed by western blot with Catalog No:117102(BCL6 Antibody) at dilution of 1:1500

IP Result of anti-BCL6 (IP:Catalog No:117102, 4ug; Detection:Catalog No:117102 1:1000) with Raji cells lysate 1400ug.

