

# BNIP3L antibody



Catalog Number: 108474

## Product name

BNIP3L antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

BNIP3L Rabbit Polyclonal antibody. Positive WB detected in mouse kidney tissue, K-562 cells. Observed molecular weight by Western-blot: 70 kDa

## Preparation

This antibody was obtained by immunization of BNIP3L recombinant protein (Accession Number: NM\_004331). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:117209(BNIP3L Antibody) at dilution of 1:600