

# CA14 antibody



Catalog Number: 108746

## Product name

CA14 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

CA14 Rabbit Polyclonal antibody. Positive WB detected in mouse eye tissue. Observed molecular weight by Western-blot: 45-50kd

## Preparation

This antibody was obtained by immunization of CA14 recombinant protein (Accession Number: NM\_012113). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



mouse eye tissue were subjected to SDS PAGE followed by western blot with Catalog No:108746(CA14 antibody) at dilution of 1:400