

# CD4 antibody

Catalog Number: 109119

## General Information

### Product name:

CD4 antibody

### Immunogen:

[Human CD4 Recombinant protein \(GST tag\)](#)

### Applications:

ELISA, IP, WB

### Specificity:

Human, Mouse, Rat; other species not tested.

### Antibody description:

CD4 Rabbit Polyclonal antibody. Positive IP detected in mouse thymus tissue. Positive WB detected in mouse thymus tissue. Observed molecular weight by Western-blot: 55 kDa

## Properties

### Preparation:

This antibody was obtained by immunization of CD4 recombinant protein (Accession Number: BC025782). Purification method: Antigen affinity purified.

### Formulation:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

### Storage:

Store at -20°C. DO NOT ALIQUOT

### Clonality:

Polyclonal

**Ig Type:** Rabbit IgG


## Applications

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

## Validations



mouse thymus tissue were subjected to SDS PAGE followed by western blot with Catalog No:109119(CD4 antibody) at dilution of 1:1000



IP Result of anti-CD4 (IP:Catalog No:109119, 4ug; Detection:Catalog No:109119 1:500) with mouse thymus tissue lysate 2800ug.