

CDON antibody



Catalog Number: 109171

Product name

Polyclonal

CDON antibody

Ig Type

Specificity

Rabbit IgG

human, mouse, rat; other species not tested.

Applications

Antibody description

ELISA, WB

CDON Rabbit Polyclonal antibody. Positive WB detected in SH-SY5Y cells, HeLa cells. Observed molecular weight by Western-blot: 140kd

Dilutions

Recommended Dilution:

Preparation

WB: 1:500-1:5000

This antibody was obtained by immunization of CDON recombinant protein (Accession Number: XM_011542863). Purification method: Antigen affinity purified.

Validations


Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:109171(CDON antibody) at dilution of 1:1000