

# CDS2 antibody



Catalog Number: 109174

## Product name

Polyclonal

CDS2 antibody

## Ig Type

## Specificity

Rabbit IgG

human, mouse, rat; other species not tested.

## Applications

## Antibody description

ELISA, WB

CDS2 Rabbit Polyclonal antibody. Positive WB detected in Y79 cells, human heart tissue. Observed molecular weight by Western-blot: 51kd

## Dilutions

Recommended Dilution:

## Preparation

WB: 1:500-1:5000

This antibody was obtained by immunization of CDS2 recombinant protein (Accession Number: NM\_003818). Purification method: Antigen affinity purified.

## Validations


## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality



Y79 cells were subjected to SDS PAGE followed by western blot with Catalog No:109174(CDS2 antibody) at dilution of 1:500