

CENPJ antibody

Catalog Number: 109254

Product name

CENPJ antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

CENPJ Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 153kd

Preparation

This antibody was obtained by immunization of CENPJ recombinant protein (Accession Number: BC024209). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:109254(CENPJ antibody) at dilution of 1:500



Immunofluorescent analysis of HepG2 cells using Catalog No:109254(CENPJ Antibody) at dilution of 1:25 and Rhodamine-Goat anti-Rabbit IgG