

# CHRNA5 antibody



Catalog Number: 109280

## Product name

CHRNA5 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

CHRNA5 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, A431 cells, A549 cells, mouse bladder tissue, mouse stomach tissue. Observed molecular weight by Western-blot: 53 kDa

## Preparation

This antibody was obtained by immunization of CHRNA5 recombinant protein (Accession Number: NM\_000745). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:109280(CHRNA5 antibody) at dilution of 1:500