

# CLEC12B antibody



Catalog Number: 109347

## Product name

CLEC12B antibody

## Clonality

Polyclonal

## Immunogen

[Human CLEC12B Recombinant protein \(His tag\)](#)

## Ig Type

Rabbit IgG

## Specificity

Human, Mouse, Rat; other species not tested.

## Applications

ELISA, WB

## Antibody description

CLEC12B Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, K-562 cells, mouse spleen tissue, rat brain tissue. Observed molecular weight by Western-blot: 32 kDa


## Dilutions

Recommended Dilution:  
WB: 1:500-1:5000

## Preparation

This antibody was obtained by immunization of CLEC12B recombinant protein (Accession Number: NM\_205852). Purification method: Antigen affinity purified.

## Validations



## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:109347(CLEC12B Antibody) at dilution of 1:1000