

Complement factor B antibody



Catalog Number: 109446

Product name

Complement factor B antibody

Specificity

Human; other species not tested.

Antibody description

Complement factor B Rabbit Polyclonal antibody. Positive IHC detected in human liver tissue. Positive WB detected in human serum tissue, human blood tissue. Observed molecular weight by Western-blot: 93-100 kDa, 60-69 kDa

Preparation

This antibody was obtained by immunization of Complement factor B recombinant protein (Accession Number: NM_001710). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



human serum tissue were subjected to SDS PAGE followed by western blot with Catalog No:109446(CFB antibody) at dilution of 1:800



Immunohistochemistry of paraffin-embedded human liver tissue slide using Catalog No:109446(CFB Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human liver tissue slide using Catalog No:109446(CFB Antibody) at dilution of 1:50

Complement factor B antibody



Catalog Number: 109446

(under 40x lens)