

CST6 antibody



Catalog Number: 109613

Product name

ELISA, IHC

CST6 antibody

Dilutions

Specificity

Recommended Dilution:

Human, Mouse, Rat; other species not tested.

IHC: 1:20-1:200


Antibody description

CST6 Rabbit Polyclonal antibody. Positive IHC detected in mouse liver.

Validations

Preparation

This antibody was obtained by immunization of CST6 recombinant protein (Accession Number: NM_001323). Purification method: Antigen affinity purified.



Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Immunohistochemical of paraffin-embedded mouse liver using Catalog No:109613(CST6 antibody) at dilution of 1:50 (under 40x lens)

Storage

Store at -20°C. DO NOT ALIQUOT


Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications



Immunohistochemical of paraffin-embedded mouse liver using Catalog No:109613(CST6 antibody) at dilution of 1:50 (under 10x lens)