

CST6 antibody

Catalog Number: 109613

Product name

CST6 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

CST6 Rabbit Polyclonal antibody. Positive IHC detected in mouse liver.

Preparation

This antibody was obtained by immunization of CST6 recombinant protein (Accession Number: NM_001323). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded mouse liver using Catalog No:109613(CST6 antibody) at dilution of 1:50 (under 40x lens)



Immunohistochemical of paraffin-embedded mouse liver using Catalog No:109613(CST6 antibody) at dilution of 1:50 (under 10x lens)