

# CXCR2 antibody



Catalog Number: 109643

## Product name

CXCR2 antibody

## Immunogen

[Human IL8RB Recombinant protein \(GST tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

CXCR2 Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive IP detected in HepG2 cells. Positive WB detected in L02 cells, COLO 320 cells. Observed molecular weight by Western-blot: 45-50 kDa

## Preparation

This antibody was obtained by immunization of CXCR2 recombinant protein (Accession Number: XM\_005246530). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF, IP

## Dilutions


Recommended Dilution:

WB: 1:1000-1:10000


IP: 1:500-1:5000

IF: 1:10-1:100

## Validations



L02 cells were subjected to SDS PAGE followed by western blot with Catalog No:109643(IL8RB antibody) at dilution of 1:2000



Immunofluorescent analysis of HepG2 cells, using IL8RB antibody Catalog No:109643 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit

# CXCR2 antibody

Catalog Number: 109643

IP Result of anti-IL8RB (IP:Catalog No:109643, 4ug; Detection:Catalog No:109643 1:1000) with HepG2 cells lysate 3200ug.

