

CTNNBL1 antibody



Catalog Number: 109726

Product name

CTNNBL1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

CTNNBL1 Rabbit Polyclonal antibody. Positive IHC detected in human testis tissue. Positive WB detected in Jurkat cells, mouse spleen tissue, mouse thymus tissue. Observed molecular weight by Western-blot: 65kd

Preparation

This antibody was obtained by immunization of CTNNBL1 recombinant protein (Accession Number: BC036739). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:109726(CTNNBL1 antibody) at dilution of



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:109726(CTNNBL1 Antibody) at dilution of 1:50 (under 40x lens)