

DAZL antibody

Catalog Number: 109874

Product name

DAZL antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

DAZL Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue. Positive IP detected in mouse testis tissue. Observed molecular weight by Western-blot: 35kd

Preparation

This antibody was obtained by immunization of DAZL recombinant protein (Accession Number: NM_001351). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:109874(DAZL antibody) at dilution of 1:400



IP Result of anti-DAZL (IP:Catalog No:109874, 4ug; Detection:Catalog No:109874 1:500) with mouse testis tissue lysate 4400ug.