

# Dermatopontin antibody



Catalog Number: 109895

## Product name

Dermatopontin antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

Dermatopontin Rabbit Polyclonal antibody.  
Positive WB detected in rat skin tissue. Observed molecular weight by Western-blot: 22 kDa

## Preparation

This antibody was obtained by immunization of Dermatopontin recombinant protein (Accession Number: NM\_001937). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



rat skin tissue were subjected to SDS PAGE followed by western blot with Catalog No:109895(DPT antibody) at dilution of 1:1000