

Dermatopontin antibody



Catalog Number: 109895

Product name

Dermatopontin antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Dermatopontin Rabbit Polyclonal antibody.
Positive WB detected in rat skin tissue. Observed
molecular weight by Western-blot: 22 kDa

Preparation

This antibody was obtained by immunization of
Dermatopontin recombinant protein (Accession
Number: NM_001937). Purification method:
Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH
7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



rat skin tissue were subjected to SDS PAGE
followed by western blot with Catalog
No:109895(DPT antibody) at dilution of 1:1000