

DNASE1L3 antibody



Catalog Number: 109959

Product name

DNASE1L3 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

DNASE1L3 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells. Observed molecular weight by Western-blot: 33kd

Preparation

This antibody was obtained by immunization of DNASE1L3 recombinant protein (Accession Number: NM_004944). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:109959(DNASE1L3 antibody) at dilution of 1:600