

DNM1L,DLP1 antibody



Catalog Number: 109962

Product name

DNM1L,DLP1 antibody

Immunogen

[Human DNM1L,DLP1 Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse; other species not tested.

Antibody description

DNM1L,DLP1 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, mouse heart tissue. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 70-78 kDa

Preparation

This antibody was obtained by immunization of DNM1L,DLP1 recombinant protein (Accession Number: NM_012063). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IP, WB


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IP: 1:200-1:2000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:109962(DNM1L,DLP1 Antibody) at dilution of



IP Result of anti-DNM1L,DLP1 (IP:Catalog No:109962, 4ug; Detection:Catalog No:109962 1:500) with mouse brain tissue lysate 4000ug.