

# DIS3L2 antibody



Catalog Number: 109977

## Product name

DIS3L2 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

DIS3L2 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, mouse embryo tissue. Observed molecular weight by Western-blot: 65 kDa-99 kda

## Preparation

This antibody was obtained by immunization of DIS3L2 recombinant protein (Accession Number: BC030113). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:109977(DIS3L2 Antibody) at dilution of 1:600