

DNAJB9 antibody

Catalog Number: 110100

General Information

Product name:

DNAJB9 antibody

Immunogen:

Applications:

ELISA, WB, IP

Specificity:

Human, Mouse, Rat; other species not tested.

Antibody description:

DNAJB9 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue. Positive IP detected in mouse testis tissue. Observed molecular weight by Western-blot: 26-30 kDa

Properties

Preparation:

This antibody was obtained by immunization of DNAJB9 recombinant protein (Accession Number: BC028912). Purification method: Antigen affinity purified.

Formulation:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage:

Store at -20°C. DO NOT ALIQUOT

Clonality:

Polyclonal

Ig Type:Rabbit IgG


Applications

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:110100(DNAJB9 Antibody) at dilution of 1:600



IP Result of anti-DNAJB9 (IP:Catalog No:110100, 4ug; Detection:Catalog No:110100 1:500) with mouse testis tissue lysate 4800ug.