

EIF3H antibody



Catalog Number: 110196

General Information

Product name:

EIF3H antibody

Immunogen:

Applications:

ELISA, WB

Specificity:

human, mouse, rat; other species not tested.

Antibody description:

EIF3H Rabbit Polyclonal antibody. Positive WB detected in NIH/3T3 cells. Observed molecular weight by Western-blot: 40kd

Properties

Preparation:

This antibody was obtained by immunization of EIF3H recombinant protein (Accession Number: BC000386). Purification method: Antigen affinity purified.

Formulation:

PBS with 0.1% sodium azide and 50% glycerol pH

7.3.

Storage:

Store at -20°C. DO NOT ALIQUOT

Clonality:

Polyclonal


Ig Type:Rabbit IgG

Applications

Recommended Dilution:

WB: 1:200-1:2000

Validations



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with Catalog No:110196(EIF3H antibody) at dilution of 1:500