

Ephrin B1 antibody



Catalog Number: 110386

Product name

Ephrin B1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Ephrin B1 Rabbit Polyclonal antibody. Positive WB detected in mouse lung tissue, human testis tissue, mouse embryo tissue. Observed molecular weight by Western-blot: 50kd

Preparation

This antibody was obtained by immunization of Ephrin B1 recombinant protein (Accession Number: NM_004429). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:110386(EFNB1 antibody) at dilution of 1:300