

# Ephrin B1 antibody



Catalog Number: 110386

## Product name

Polyclonal

Ephrin B1 antibody

## Ig Type

## Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

## Applications

## Antibody description

ELISA, WB

Ephrin B1 Rabbit Polyclonal antibody. Positive WB detected in mouse lung tissue, human testis tissue, mouse embryo tissue. Observed molecular weight by Western-blot: 50kd

## Dilutions


Recommended Dilution:

WB: 1:200-1:2000

## Preparation

This antibody was obtained by immunization of Ephrin B1 recombinant protein (Accession Number: NM\_004429). Purification method: Antigen affinity purified.

## Validations



## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:110386(EFNB1 antibody) at dilution of 1:300