

# EZH1 antibody



Catalog Number: 110424

## Product name

EZH1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

EZH1 Rabbit Polyclonal antibody. Positive WB detected in RAW 264.7 cells, HeLa cells, HL-60 cells, RAW264.7 cells. Observed molecular weight by Western-blot: 77-86 kDa, 69 kDa

## Preparation

This antibody was obtained by immunization of EZH1 recombinant protein (Accession Number: NM\_001991). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

## Validations



RAW264.7 cells were subjected to SDS PAGE followed by western blot with Catalog No:110424(EZH1 antibody) at dilution of 1:300