

ETFB antibody



Catalog Number: 110484

Product name

ETFB antibody

Immunogen

[Human ETFB Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ETFB Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive WB detected in mouse skeletal muscle tissue, human heart tissue, mouse brain tissue, rat brain tissue. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 28 kDa

Preparation

This antibody was obtained by immunization of ETFB recombinant protein (Accession Number: NM_001985). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IF, WB, IP

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


IP: 1:500-1:5000

IF: 1:10-1:100

Validations



Immunofluorescent analysis of HepG2 cells, using ETFB antibody Catalog No:110484 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:110484(ETFB Antibody) at dilution of 1:1000

ETFB antibody

Catalog Number: 110484

IP Result of anti-ETFB (IP:Catalog No:110484, 3ug; Detection:Catalog No:110484 1:1000) with mouse brain tissue lysate 4000ug.

