

# ETFB antibody



Catalog Number: 110484

## Product name

ETFB antibody

IP: 1:500-1:5000

IF: 1:10-1:100

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

ETFB Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive WB detected in mouse skeletal muscle tissue, human heart tissue, mouse brain tissue, rat brain tissue. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 28 kDa

## Preparation

This antibody was obtained by immunization of ETFB recombinant protein (Accession Number: NM\_001985). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IF, WB, IP

## Dilutions


Recommended Dilution:

WB: 1:500-1:5000


## Validations



Immunofluorescent analysis of HepG2 cells, using ETFB antibody Catalog No:110484 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:110484(ETFB Antibody) at dilution of 1:1000



IP Result of anti-ETFB (IP:Catalog No:110484, 3ug; Detection:Catalog No:110484 1:1000) with

# ETFB antibody



Catalog Number: 110484

---

mouse brain tissue lysate 4000ug.