

FANCL antibody



Catalog Number: 110521

Product name

FANCL antibody

Immunogen

[Human FANCL Recombinant protein \(His tag\)](#)

Specificity

human, mouse, rat; other species not tested.

Antibody description

FANCL Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, HeLa cells. Observed molecular weight by Western-blot: 41-43kd

Preparation

This antibody was obtained by immunization of FANCL recombinant protein (Accession Number: NM_018062). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:2000

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:110521(FANCL antibody) at dilution of 1:300