

# FECH antibody



Catalog Number: 110613

---

## Product name

FECH antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

FECH Rabbit Polyclonal antibody. Positive WB detected in human heart tissue, HeLa cells, human liver tissue, K-562 cells. Observed molecular weight by Western-blot: 46kd

## Preparation

This antibody was obtained by immunization of FECH recombinant protein (Accession Number: NM\_000140). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:110613(FECH antibody) at dilution of 1:500