

# FKBP7 antibody



Catalog Number: 110680

## Product name

ELISA, WB, IHC

FKBP7 antibody

## Dilutions

## Specificity

Recommended Dilution:

Human, Mouse, Rat; other species not tested.

WB: 1:500-1:5000

IHC: 1:20-1:200


## Antibody description

FKBP7 Rabbit Polyclonal antibody. Positive IHC detected in human lung cancer tissue. Positive WB detected in human lung tissue, human heart tissue. Observed molecular weight by Western-blot: 30 kDa

## Validations

## Preparation

This antibody was obtained by immunization of FKBP7 recombinant protein (Accession Number: NM\_181342). Purification method: Antigen affinity purified.



## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

human lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:110680(FKBP7 antibody) at dilution of 1:500

## Storage

Store at -20°C. DO NOT ALIQUOT


## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications



Immunohistochemical of paraffin-embedded human lung cancer using Catalog No:110680(FKBP7 antibody) at dilution of 1:50 (under 10x lens)