

FKBP7 antibody



Catalog Number: 110680

Product name

FKBP7 antibody

Immunogen

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

FKBP7 Rabbit Polyclonal antibody. Positive IHC detected in human lung cancer tissue. Positive WB detected in human lung tissue, human heart tissue. Observed molecular weight by Western-blot: 30 kDa

Preparation

This antibody was obtained by immunization of FKBP7 recombinant protein (Accession Number: NM_181342). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



human lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:110680(FKBP7 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human lung cancer using Catalog No:110680(FKBP7 antibody) at dilution of 1:50 (under 10x lens)